



# TRILUTION® micro Software v3.0

## INSTALLATION QUALIFICATION PROCEDURES

The installation qualification (IQ) procedures are intended for end users who are implementing good laboratory practice (GLP) requirements and can be readily incorporated into proprietary standard operating procedures (SOPs).

### 1. Product Identification

#### Manufacturer

Gilson, Inc.  
3000 Parmenter Street  
PO Box 620027  
Middleton, WI 53562-0027 USA

Name:	TRILUTION® micro Software
Version:	3.0

#### Supplier

Organization	<hr/>
Address	<hr/> <hr/> <hr/>
Phone Number	<hr/>
Fax Number	<hr/>
Name of Gilson Representative	<hr/>

## User

Organization	<hr/>
Department	<hr/>
Site (Room)	<hr/>
Date of Installation	<hr/>
System that Includes the Software	<hr/>
<b>Primary Contact</b>	
Name	<hr/>
Phone	<hr/>
Email	<hr/>

## 2. Installation

TRILUTION® micro Software is factory-installed on the tablet.


Follow the instructions provided in *TRILUTION micro v3.0 Installation Instructions* (part number LT375091) that describes how to install TRILUTION® micro Software v3.0 on a PC.

### Documents

In addition to this document, the following are provided. For reference purposes, specify the storage location for each.

TRILUTION® micro Software Documentation USB	<hr/>
TRILUTION® micro Overview	<hr/>
TRILUTION® micro Installation Instructions	<hr/>

## TRILUTION® micro Software

1. Start TRILUTION micro. Select **Start** and then select **(All apps) > Gilson Applications > TRILUTION micro 3.0**.
2. Select  to access the **About** screen that lists the software version and then select **System Info**.
3. Optionally, select **Export** to save the system info to a file or **Clipboard** to paste it into a file or e-mail.

Name	Observed Version	Expected Version
Gilson TRILUTION micro	_____	3.0.18.1
Gilson TRILUTION micro Execution Engine	_____	1.3.13.0
Gilson TRILUTION micro Embedded Application	_____	1.4.75.0
Gilson TRILUTION micro Labware Creator	_____	1.4.8.0
Gilson TRILUTION micro GBS Interface	_____	1.4.2.0
Gilson TRILUTION micro Normalization Assistant	_____	1.4.2.0
Gilson TRILUTION micro qPCR Assistant	_____	1.4.2.0

## PIPETMAX Alignment Utility

1. Close TRILUTION micro. The PIPETMAX Alignment Utility will not start if TRILUTION micro is running.
2. Click the **Start** button and then select: **(All apps) > Gilson Applications > PIPETMAX Alignment Utility**.

The version is displayed in the lower, left corner of the PIPETMAX Alignment Utility.

Name	Observed Version	Expected Version
PIPETMAX Alignment Utility	_____	1.0.77.0

### 3. Acceptance

If no installation problems were detected or installation problems that were detected were resolved, have the local Gilson representative who verified the installation of TRILUTION micro provide the information requested below.

Gilson Representative	_____
Organization	_____
Signature	_____
Date	_____
<input type="checkbox"/> End user was provided with a copy of this document.	