




Purification of a 25-mer ssDNA Oligonucleotide using the PLC 2020 Personal Purification System



Conditions

Sample: 25-mer ssDNA oligonucleotide, DMT-On (single-stranded, Dimethoxytrityl-On), with fluorescence & non-fluorescence tags

Load: 1 optical density

HPLC column: 250 x 4.6 mm C4

Mobile phase:

- A: 0.1 M total essential amino acid pH 7
Triethyl ammonium acetate
- B: Acetonitrile

Wavelength: 260 & 550 nm

A C4 reverse phase column was used to separate a 25-mer ssDNA oligonucleotide using the Gilson PLC 2020 Personal Purification System prior to DMT-On purification. The oligonucleotide eluted at 18.038 minutes. Final purification was performed using volume collection (0.5 mL fractions). Purification was triggered on absorption at 260 nm.

To discover more Application Briefs and notes from Gilson, visit www.gilson.com/en/AI/Applications/