



# Gilson Quality Control Report

## Pipette Head Addendum

### Pipette Head for PIPETMAX®

#### Default Calibration Values

##### MAX8X20

Volume Setting	1 µL Minimum	10 µL 50%	20 µL Nominal
Calibration Volume Adjustment	1000	10000	20000
Calibration Offset	990	8140	16140

##### MAX8X200

Volume Setting	20 µL Minimum	100 µL 50%	200 µL Nominal
Calibration Volume Adjustment	20000	100000	200000
Calibration Offset	1660	8040	15960

##### MAX1000

Volume Setting	100 µL Minimum	500 µL 50%	1000 µL Nominal
Calibration Volume Adjustment	99500	496000	994000
Calibration Offset	1680	8080	16064

##### MAX4X1200

Volume Setting	120 µL Minimum	600 µL 50%	1200 µL Nominal
Calibration Volume Adjustment	120000	600000	1200000
Calibration Offset	1982	9669	19264



Refer to the *PIPETMAX® 268 User's Guide* (part number LT255519), **Add Pipette Head(s)** section for information about scanning or entering pipette head calibration values.