



Tubing Selection Guide

MINIPULS® 3 Peristaltic Pump

Tubing for MF1 and MF4 Pump Heads Length = Approximately 40 cm (16 in.)							Tubing for HF2 Pump Head Length = 3 m (118 in.)					
Material	Internal Diameter (mm)	Reference Tubing (Calibrated Tubing)	Color Coded Stops		Extension Tubing 3 m (118 in.)	Connectors	Maximum Flow Rate (mL/min)	Material	Internal Diameter (mm)	Reference Tubing (Calibrated Tubing)	Maximum Flow Rate (mL/min)	
PVC	0.25	F117932	●	●	F117952	F117985	0.574	PVC	2.06	F117965	33	
	0.38	F117933	●	●	F117953	F117986	1.22		2.29	F117966	41	
	0.51	F117934	●	●	F117954	F117986	1.90		2.79	F117968	66	
	0.76	F117936	●	●	F117956	F117987	3.00		3.16	F117969	85	
	1.02	F117938	○	○	F117958	F1179941	6.0		4.0	F117970	128	
	1.30	F117940	○	○	F117960	F1179941	9.0		5.0	F117980	179	
	1.52	F117942	●	●	F117962	F1179941	14		6.0	F117981	218	
	1.65	F117943	●	●	F117963	F1179941	16		8.0	F117982	340	
	2.06	F117945	●	●	F117965	F1179941	25		Iso-Versinic® or Viton®	2.0	F117744	30
	2.29	F117946	●	●	F117966	F1179941	29			3.0	F117746	67
	2.79	F117948	●	○	F117968	F1179951	39			4.0	F117748	120
	3.18	F117949	●	○	F117969	F1179951	47			6.0	F117750	246
Iso-Versinic® or Viton®	0.5	F1817741	●	●	F117740	F117986	2.9	Silicone	2.0	F1825133	31	
	1.0	F1817743	●	●	F117742	F1179941	10		3.0	F117975	57	
	2.0	F1817745	●	●	F117744	F1179941	24		4.0	F117976	105	
	3.0	F117747	●	●	F117746	F1179951	55		5.0	F117977	148	
	4	F117749	●	●	F117748	-	70		6.0	F117978	207	
Silicone	0.64	F1825111	●	○	F1825131	F117987	3.2	7.0	F117979	300		
	1.02	F1825112	○	○	F1825132	F1179941	7	Santoprene™ or PharMed®	2.0	F1825123	31	
	2.06	F1825113	●	●	F1825133	F1179951	21		3.2	F1825124	77	
	2.79	F1825114	●	○	F117975	F1179951	31		6.4	F1825125	240	
Santoprene™ or PharMed®	0.51	F1825101	●	●	F1825121	F117986	2.5					
	1.02	F1825102	○	○	F1825122	F1179941	7					
	2.06	F1825103	●	●	F1825123	F1179941	24					

All product and company names are trademarks™ or registered® trademarks of their respective holders. Use of the trademark(s) in this document does not imply any affiliation with or endorsements by the trademark holder(s).