



pipetman® 

Volume and Increment Table

<i>PIPETMAN® M model</i>	<i>Part Number</i>	<i>Volume Range Pipet Mode</i>	<i>Volume Range Repetitive Mode</i>	<i>Maximum Number of Aliquots (Repetitive Mode)</i>	<i>Increment</i>
P10M	F81022	0.5–10 µL	0.5–10 µL	20	0.01 µL
P20M	F81023	2–20 µL	2–20 µL	10	0.01 µL
P100M	F81013	5–100 µL	5–100 µL	20	0.1 µL
P200M	F81024	20–200 µL	5–200 µL	40	0.1 µL
P300M	F81014	20–300 µL	10–300 µL	30	0.1 µL
P1200M	F81015	100–1200 µL	20–1200 µL	60	1 µL
P5000M	F81016	500–5000 µL	100–5000 µL	50	1 µL
P10mL M	F81017	1–10 mL	200 µL–10 mL	50	0.01 mL
P8x10M	F81025	0.5–10 µL	0.5–10 µL	20	0.01 µL
P12x10M	F81026	0.5–10 µL	0.5–10 µL	20	0.01 µL
P8x20M	F81027	1–20 µL	1–20 µL	20	0.01 µL
P12x20M	F81028	1–20 µL	1–20 µL	20	0.01 µL
P8x100M	F81018	10–100 µL	5–100 µL	20	0.1 µL
P12x100M	F81019	10–100 µL	5–100 µL	20	0.1 µL
P8x200M	F81029	20–200 µL	5–200 µL	40	0.1 µL
P12x200M	F81030	20–200 µL	5–200 µL	40	0.1 µL
P8x300M	F81031	10–300 µL	10–300 µL	30	0.1 µL
P12x300M	F81032	10–300 µL	10–300 µL	30	0.1 µL
P8x1200M	F81020	50–1200 µL	50–1200 µL	24	1 µL
P12x1200M	F81021	50–1200 µL	50–1200 µL	24	1 µL

Should you have any questions regarding PIPETMAN M, please contact Nora Meneceur, MLH Product Manager, at nmeneceur@gilson.com.



Gilson, Inc. | PO Box 620027 | Middleton, WI 53562-0027 | Tel: 608-836-1551 OR 800-445-7661 | Fax: 608-831-4451

Gilson S.A.S. | 19, avenue des Entrepreneurs | BP 145, F-95400 Villiers-le-bel, FRANCE

www.gilson.com | sales@gilson.com | service@gilson.com | training@gilson.com

LT801574-01