


Guaranteed Performance with Outstanding Accuracy
ELECTRONIC PIPETTES



| On-Target Performance

Trust in your results with guaranteed* performance in both standard and repetitive pipetting modes.

| Pipetting Comfort

Lightweight and perfectly balanced with low ejection and volume setting forces, for hours of pipetting comfort.

| Ease-of-Use

With a streamlined user interface and five pipetting modes, control of your sample prep is at your fingertips. Switch to custom mode to create a new protocol. Save and transfer your protocol securely into PIPETMAN M software for recall at a later date.

*For more information about guaranteed performance, please visit www.gilson.com



PIPETMAN M Single Channel

Model	Part Number	Volume Range		PIPETMAN* DIAMOND Tips	
		STANDARD PIPET Mode	REPETITIVE PIPET Mode		
P10M	F81022	0.5 - 10 µL	0.5 - 10 µL	D10* DL10*	DF10ST DFL10ST
P20M	F81023	2 - 20 µL	2 - 20 µL	D200	DF30ST
P100M	F81013	5 - 100 µL	5 - 100 µL	D200	DF100ST
P200M	F81024	20 - 200 µL	5 - 200 µL	D200 D300	DF200ST DF300ST
P300M	F81014	20 - 300 µL	10 - 300 µL	D200** D300	DF200ST** DF300ST
P1200M	F81015	100 - 1200 µL	20 - 1200 µL	D1000** D1200	DF1000ST** DF1200ST
P5000M	F81016	500 - 5000 µL	100 - 5000 µL	D5000	
P10mLM	F81017	1 - 10 mL	200 µL - 10 mL	D10mL	

PIPETMAN M Multichannel

P8x10M	F81025	0.5 - 10 µL	0.5 - 10 µL	D10 DL10	DF10ST DFL10ST
P12x10M	F81026				
P8x20M	F81027	1 - 20 µL	1 - 20 µL	DL10 D200	DF30ST
P12x20M	F81028				
P8x100M	F81018	10 - 100 µL	5 - 100 µL	D200	DF100ST
P12x100M	F81019				
P8x200M	F81029	20 - 200 µL	5 - 200 µL	D200 D300	DF100ST** DF200ST DF300ST
P12x200M	F81030				
P8x300M	F81031	10 - 300 µL	10 - 300 µL	D200** D300	DF200ST** DF300ST
P12x300M	F81032				
P8x1200M	F81020	50 - 1200 µL	50 - 1200 µL	D1200	DF1200ST
P12x1200M	F81021				



EMPOWER YOUR PIPETMAN M

Accessories	Part Number
PIPETMAN M Power Carrousel	FB1001
PIPETMAN M Stand Adaptor	F807023

Gilson maximum permissible errors are guaranteed only when PIPETMAN pipettes are used with the recommended PIPETMAN DIAMOND Tips. Visit www.gilson.com

Under these conditions, Gilson volumetric specifications in standard pipetting (PIPET mode) are guaranteed with a performance exceeding ISO 8655-2 recommendations for this mode.

* A plastic adapter is required to eject D10 and DF10ST tips (short tips). The adapter is supplied with P10M. No adapter is required when using DL10 and DFL10ST tips (long tips).

** These PIPETMAN DIAMOND Tips can be used with the indicated pipettes for the REPETITIVE mode but not until the maximal volume of the pipette in PIPET mode. Please refer to the volume range of your tips.

Recommended tips to get best results in REPETITIVE mode.

www.gilson.com/contactus

ISO 9001
CERTIFIED

LT800593/A | ©2019 Gilson S.A.S.
All rights reserved. All trademarks are the property of Gilson, Inc. and its subsidiaries.

 **GILSON**[®]